

Centrol, Inc.

24/7 Medical Alert Personal Emergency Response System **INSTALLATION AND TESTING PROCEDURE**

Installation:

1. Select a convenient location that is next to both a telephone jack and a household electrical outlet. The kitchen is not a good location because of proximity to appliances, which may generate interference.
2. Plug the power adapter into an active electrical outlet. Once plugged in, the Unit will “beep” alerting you that the system is now activated and the power indicator light will light. CAUTION: Be sure this outlet can not be turned off by a wall switch.
3. If there is a phone plugged in the selected phone jack, unplug it. Then plug the telephone line cord from the PERS Unit labeled “THIS JACK TO WALL OUTLET” into the phone jack. Now you may plug your phone or other device into the jack on the back of the PERS Unit labeled “PLUG PHONE OR OTHER DEVICE HERE.”

Battery Information:

1. Your PERS Unit comes with a rechargeable back-up battery inside the housing. The Unit must be plugged into your electrical outlet for up to 24 hours to assure full back-up battery operation. If your electrical power goes out, your PERS will remain operational for approximately 8+ hours, powered by the factory installed maintenance free rechargeable battery. Do not remove the power cord from the outlet while in use.
2. Pendant transmitters’ lithium battery has an expected life of up to 10 years or 50,000 activations, whichever comes first.

Testing:

It is recommended to test your unit upon installation **AND EVERY 30 DAYS**. Take the transmitter (pendant) in hand and walk away from the PERS Unit. Press the transmitter. The PERS will sound (acknowledging it received the signal from the transmitter). Let the call go through and when you talk with the professional, simply say... “I am just testing. There is no emergency. Thank you.” The test is finished.

Remember, Your transmitter will work only in and around your residence. And, if you move, or change phone numbers, call CENTROL with your new information BEFORE you move.

IMPORTANT: If you ship or store your system, you will need to Unplug the power supply and press and hold the RESET button. When you hear two (2) beeps, release the reset button. The backup battery will then be disconnected. The Unit will then be safe to ship. If you do not power down, the unit will begin to beep telling you the back-up battery is about to die. Please be sure to take this action especially when shipping either UPS or US Postal Service.

Centrol, Inc. Emergency Response Systems
410-828-5564 or 877-792-4494 toll free
FAX: 410-828-7730